LOT DBCL18041
DBCL18042

 2018-04-30

Thermo
SCIENTIFIC

MAS[®] Diabetes

LIQUID ASSAYED DIABETES CONTROL

IVD

INTENDED USE

MAS[®] Diabetes is intended for use as an assayed control material for monitoring Glycosylated Hemoglobin (A1c) assay procedures.

PRODUCT DESCRIPTION

Diabetes is an assayed product prepared from human whole blood adjusted to specific concentrations of glycated hemoglobin. Preservatives and stabilizers are added to maintain product integrity.

CAUTION: Diabetes is prepared from human source material. Components of the control which are derived from human source material have been tested using FDA accepted methods and found non-reactive for Hepatitis B Surface Antigen (HBsAg), Hepatitis C (HCV), HIV-1 and HIV-2. However, no test method can offer complete assurance that products derived from human source material are free of infectious agents. This control must be handled in accordance with recommendations from Centers for Disease Control/National Institutes of Health manual, "Biosafety in Microbiological and Biomedical Laboratories," 2009.

CAUTION: The packaging of this product contains dry natural rubber.

STORAGE AND STABILITY

Diabetes stored at -25 to -15°C is stable until the expiration date on the box. Once opened, vials of **Diabetes** are stable for 30 days when stored tightly capped at 2-8°C. Unopened vials of **Diabetes** are stable for 60 days from receipt when stored at 2-8°C. The control can be thawed and refrozen one time. **Self-defrosting freezers are not suitable.**

Bacterial contamination produces an increase in turbidity and/or a characteristic odor. Discard vial if evidence of microbial contamination is observed.

CONTROL RANGES

The assigned ranges for these controls are based upon replicate assays of representative samples of the product by participating laboratories in accordance with established protocol. All values have been assigned with instruments and instrument manufacturer's reagents available at the time of assay. Subsequent instrument or reagent modifications may invalidate these assigned ranges.

Expected values may vary slightly with different reagent and/or methodologies used. Refer to the included table for values obtained for specific systems. Values listed are specific for this lot of control only. Good laboratory practice suggests that each laboratory establish its own parameters.

INSTRUCTIONS FOR USE

Thaw control at room temperature (18-25°C) on a rocker or with periodic gentle inversion until liquid and then immediately store at 2-8°C. Thoroughly mix the contents of the vial before each use by gently inverting for several minutes. Open the vial and transfer the required quantity of control into a clean sample cup. Replace cap immediately and store the opened vial at 2-8°C. Assay controls in accordance with the reagent manufacturer's assay's and/or instrument's instructions for unknown specimens.

QUALITY CONTROL

All quality control requirements should be performed in conformance with local, state and/or federal regulations or accreditation requirements.

LIMITATIONS OF PROCEDURE

Values for each lot are given in the product insert. Different laboratories, however, may expect slightly different results. It is recommended that each laboratory establish its own expected range.

Accurate and reproducible results are dependent upon properly functioning instruments, reagents, and good laboratory technique. This product is intended for use as an assayed control for quantitative assays of listed constituents in human whole blood. This product is not intended for use as a calibrator. For professional use only.

TECHNICAL ASSISTANCE

In the USA, for technical assistance, call 800-232-3342 or 510-979-5417. For insert updates and information, if your laboratory subscribes to LabLink xL, visit www.maslablink.com and select LabLink Extra. Alternatively, to subscribe to LabLink xL call 800-232-3342 or 510-979-5451.

Outside of the USA, if your laboratory subscribes to LabLink xL, visit www.maslablink.com select LabLink Extra. Alternatively, please contact your local sales office or authorized distributor.

Cat. No.
DBCL-MP

Description
MAS[®] Diabetes, Multi-Pack

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USA
☎ 800-232-3342
sales.diagnostics.fmt@thermofisher.com

China
☎ +86 800-810-5118
cdx.cn.info@thermofisher.com

New Zealand
☎ 0800 933 966
auinfo@thermofisher.com

Switzerland & Austria
☎ +41 26 663 86 70
cdx.ch.info@thermofisher.com

Asia Pacific
☎ +61 1800 333 110
cdd.asia.info@thermofisher.com

France
☎ +33 1 40 86 65 20
cdx.fr.info@thermofisher.com

Nordic
☎ +47 2 325 0433
info.nordic.cdd@thermofisher.com

United Kingdom & Ireland
☎ +44 1442 868 940
cdx.uk.info@thermofisher.com

Australia
☎ +61 1800 333 110
auinfo@thermofisher.com

Germany
☎ +49 0800-40 40 771
cdx.de.info@thermofisher.com

South Africa
☎ +27117926790
support-za.idd@thermofisher.com

Canada
☎ 800-282-4075
info.cddcanada@thermofisher.com

Japan
☎ +81 (0)120-147-075
JPYOK-CDD.QC@thermofisher.com

Spain, Portugal & Italy
☎ +34 93589 8338
cdx.es.info@thermofisher.com

For countries not listed:
distributor.cdd@thermofisher.com

Microgenics Corporation
46500 Kato Road
Fremont, CA 94538 USA
ISO 13485 Certified Company
www.thermoscientific.com/diagnostics

EC REP

B.R.A.H.M.S GmbH
Neuendorfstrasse 25
16761 Hennigsdorf, Germany
Tel: +49 (0) 800 404 077 12
Fax: +49 (0) 800 404 077 13

DBCL-INS-VA
Rev. 14 2015 11

Thermo
SCIENTIFIC